

Tightening torques ZERAMEX® VICARBO® screws



System ZERAMEX® P6

   <p>P15450</p> <p>Tightening torque: 25Ncm</p>	  <p>P15451</p> <p>Tightening torque: 25Ncm</p>	   <p>P35550</p> <p>Tightening torque: 15Ncm</p>
   <p>P16450</p> <p>Tightening torque: 35Ncm</p>	  <p>P16451</p> <p>Tightening torque: 35Ncm</p>	   <p>P36550</p> <p>Tightening torque: 15Ncm</p>

System ZERAMEX® XT

     <p>RB16550</p> <p>Tightening torque: 25Ncm</p>	    <p>RB36553</p> <p>Tightening torque: 5Ncm</p>	   <p>RB36550</p> <p>Tightening torque: 15Ncm</p>
--	---	--